

AD_____

Award Number: DAMD17-98-1-8054

TITLE: Dynamic Functional Mammoscopy: A Non-Ionizing Imaging
Technique Enhancing Early Detection of Breast Cancer

PRINCIPAL INVESTIGATOR: Suzanne Smith

CONTRACTING ORGANIZATION: Columbia University
New York, New York 10032

REPORT DATE: October 2000

TYPE OF REPORT: Annual

PREPARED FOR: U.S. Army Medical Research and Materiel Command
Fort Detrick, Maryland 21702-5012

DISTRIBUTION STATEMENT: Approved for Public Release;
Distribution Unlimited

The views, opinions and/or findings contained in this report are those of the author(s) and should not be construed as an official Department of the Army position, policy or decision unless so designated by other documentation.

20020401 032

REPORT DOCUMENTATION PAGE			Form Approved OMB No. 074-0188	
Public reporting burden for this collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing this collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden to Washington Headquarters Services, Directorate for Information Operations and Reports, 1215 Jefferson Davis Highway, Suite 1204, Arlington, VA 22202-4302, and to the Office of Management and Budget, Paperwork Reduction Project (0704-0188), Washington, DC 20503				
1. AGENCY USE ONLY (Leave blank)	2. REPORT DATE October 2000	3. REPORT TYPE AND DATES COVERED Annual (1 Oct 99 - 30 Sep 00)		
4. TITLE AND SUBTITLE Dynamic Functional Mammoscopy: A Non-Ionizing Imaging Technique Enhancing Early Detection of Breast Cancer		5. FUNDING NUMBERS DAMD17-98-1-8054		
6. AUTHOR(S) Suzanne Smith				
7. PERFORMING ORGANIZATION NAME(S) AND ADDRESS(ES) Columbia University New York, New York 10032 E-Mail: sjs9@columbia.edu		8. PERFORMING ORGANIZATION REPORT NUMBER		
9. SPONSORING / MONITORING AGENCY NAME(S) AND ADDRESS(ES) U.S. Army Medical Research and Materiel Command Fort Detrick, Maryland 21702-5012		10. SPONSORING / MONITORING AGENCY REPORT NUMBER		
11. SUPPLEMENTARY NOTES				
12a. DISTRIBUTION / AVAILABILITY STATEMENT Approved for Public Release; Distribution Unlimited.			12b. DISTRIBUTION CODE	
13. ABSTRACT (Maximum 200 Words) The imaging technology of Dynamic Functional Optical Mammoscopy (DFOM) has been used to scan patients scheduled for biopsy of breast lesions. These patients were scheduled for core or excisional breast biopsy on the basis of equivocal mammographic and ancillary clinical findings within ACR BI-RADST TM categories 3 or 4. Analysis of test results of 47 patients shows that the DFOM detected cancer in nine of the 11 patients in whom biopsies confirmed malignant lesions giving a sensitivity of 82%. DFOM also correctly identified 24 of 36 benign lesions giving a specificity of 67%. In clinical practice, the adjunctive use of DFOM would have decreased the percentage of biopsies that turn out to be benign from 36/47 (77%) to 12/47 (26%). The negative predictive value, the chance that a negative DFOM result truly indicates a benign lesion, was 24/26 (92%).				
14. SUBJECT TERMS Breast Cancer			15. NUMBER OF PAGES 19	
			16. PRICE CODE	
17. SECURITY CLASSIFICATION OF REPORT Unclassified	18. SECURITY CLASSIFICATION OF THIS PAGE Unclassified	19. SECURITY CLASSIFICATION OF ABSTRACT Unclassified	20. LIMITATION OF ABSTRACT Unlimited	

Table of Contents

Front Cover	1
Report Documentation Page	2
Introduction	4
Body	4
Key Research Accomplishments	7
Reportable Outcomes	7
Conclusions	8
References	8
Appendices	9

Introduction

The imaging technology of Dynamic Functional Optical Mammoscopy (DFOM) is a breast scan based upon transmission/absorption of infrared light, which measures the dynamic patterns of breast reactivity of various physiological states in response to soft pressure. DFOM produces a functional rather than a morphological image and the dynamic pattern of tissue reactivity after mild compression is charted. Pilot study results suggested that this innovative DFOM imaging technique has the potential to determine which of the mammographically and clinically indeterminate lesions are benign vs. carcinoma and distinguish those lesions thereby avoiding biopsy. The purpose of the study reported here is to extend the preliminary results of the pilot study at Columbia Presbyterian Medical Center, using DFOM between mammography and biopsy to further evaluate the efficacy of DFOM in evaluation of breast lesions, using biopsy results to confirm diagnosis.

Body

A total of 117 patients scheduled for biopsy were scanned with the identical protocol between June 1, 2000 and September 30, 2000. The study was performed on women scheduled for core or excisional breast biopsy on the basis of equivocal mammographic and ancillary clinical findings within ACR BI-RADS™ categories 3 or 4. Women who met the selection criteria were enrolled from the normal caseload from both screening and diagnostic mammography. Each woman signed an informed consent prior to being scanned.

The scan procedure required approximately 5 minutes. During examination, the breast was placed in the soft breast holder of the system. The breast was then softly compressed by a thin transparent silicone rubber membrane using an applied pressure of approximately 10 mm Hg. For each scan, the breast was symmetrically centered on the illuminator. When the breast was correctly positioned, illumination adjustment and image recording took place following the requirements of the pressure profile.

Optical illumination was provided by an array of red light emitting diodes (LEDs) attached to the bottom surface of the soft breast holder. Light transmitted through the breast was recorded as a temporal sequence for approximately 30 seconds by a highly sensitive digital CCD camera. The image sequences were accumulated in digital memory and processed by proprietary software to accentuate differences in the temporal variations of intensity between normal/benign and malignant tissue.

Each woman was scanned by a trained technologist prior to biopsy. The scans were read by an experienced reader trained in interpreting the scans. Results were reported as either a recommendation for biopsy or a recommendation that the woman be sent to interval follow-up.

Recommendations on the basis of DOFM were compared to pathology reports of malignant or benign which were used as the gold standard. Sensitivity, specificity, and negative predictive value were calculated.

Table 1 gives the patient accounting for the 117 patients reported on during this period.

TABLE 1 – Patients Scanned June 1 – September 30, 2000

Excluded patients*	= 40
Unacceptable scans**	= 20
Scans to be interpreted	= 10
Interpreted scans	= 47
Total scans performed	= 117

*Patients who did not meet the selection criteria for the development study protocol.

**Patients whose scans were not acceptable to be interpreted.

Table 2 presents the reasons for scans determined to be unacceptable.

<u>TABLE 2 – Unacceptable Scans</u>	
Lesion sub areolar or located where it cannot be properly illuminated	= 5
Inadequate illumination in area of pathology	= 6
Device related	= 4
Incorrect illumination	= 2
Not enough breast tissue in holder	= 2
Excessive patient movement	= 1
Total	= 20

Below, Table 3 lists reasons for excluding from our study selected patients.

<u>TABLE 3 – Excluded Patients</u>	
Lesion sub areolar or located where it cannot be properly illuminated	= 12
Previous surgery in ipsilateral breast	= 9
BI-RADs 5 lesion	= 4
Post-menopausal with palpable lesion	= 3
Relevant record(s) not available at site of review	= 2
Small breast that could not be properly positioned	= 2
Patient could not remain still	= 1
Biopsy of ipsilateral breast within 3 months	= 1
Biopsy to be performed at another facility	= 1
Not recorded	= 5
Total	= 40

Table 4 presents patient demographics of race and age.

TABLE 4 – Patient Demographics

Race	
White	= 58
African American	= 5
Hispanic White	= 34
Hispanic Black	= 16
Asian	= 2
Other	= 2
Total	= 117

AGE: Average = 56
Range = 35-82

Table 5 presents the results of scan interpretation.

TABLE 5 – Results of 47 Patients Scheduled for Biopsy

		Pathology		
		Malignant	Benign	Total
DFOM Recommendation	Biopsy	9	12	21
	Interval Follow-Up	2	24	26
	Total	11	36	47

Sensitivity: 9/11 (82%)
Specificity: 24/36 (67%)
Negative Predictive Value (NPV): 24/26 (92%)

The analysis of test results on the 47 patients with interpreted scans shows that the DFOM detected cancer in 9 of the 11 patients in whom biopsies confirmed malignant lesions ("true positives"). This results in a sensitivity of 82%. The system also correctly identified 24 of 36 benign lesions ("true negatives"). In other words, the specificity of the DFOM is 24/36 (67%).

Key Research Accomplishments

Below we summarize our accomplishments in the context of the original Statement of Work, as described in the proposal.

TASK 1. Develop substructure for implementing study (months 1-3).

- (a) Obtained Institutional Review Board (IRB) Approval, including patient consent form.
- (b) Trained three consecutive research associates in protocol and procedures:
Patricia Ogiliva 11/30/98 – 6/3/99,
Homyra Hadavand 6/14/99 – 11/30/99 and
Behnaz Mesbah 12/13/99 – 8/31/01.
- (c) Develop database program to enable data entry and to allow compilation of results.

TASK 2. Develop a database of DFOM dynamic images and signatures of various pathologic lesions (months 3-34).

- (a) Recruited patients undergoing evaluation of lesions to have DFOM scan prior to tissue sampling or biopsy.
- (b) Correlated DFOM information evaluated for acquired dynamic images and signature types with pathologic lesion by type, grade (Bi-Rads scale) and size. Compared appearance with parenchymal density.

TASK 4. Interim analyses (months 12-26).

- (a) Interim analysis has been presented in this report for the first phase of the study, and is currently being collected for analysis in the next annual report.

Reportable Outcomes

The effectiveness of the DFOM in discriminating between benign and malignant breast lesions has been further evaluated and the preliminary results have been corroborated.

Conclusions

While the number of patients reported on is small, the indications of effectiveness are very encouraging. The results reported above indicate that the adjunctive use of the DOFM in clinical practice would have decreased the percentage of biopsies that turn out to be benign from 36/47 (77%) to 12/47 (26%). The negative predictive value, the chance that a negative DFOM result truly indicates a benign lesion, is 24/26 (92%).

Thus, we are encouraged by these early findings. We will continue to accrue patients and analysis data during the next phase.

References

None.

Appendices

Summary Tables of Scanned Patients - First Year of Study

CPMC

	Patient ID	Protocol	Press. P	Acceptable	Blind	Adjunct	BiRads	Unblind	Biopsy	S.conf
1	19980915A-L	post-biopsy	E	NO - NC						y
2	19980915A-R	post-biopsy	E	NO - NC						y
3	19981001A R	8408	E	YES			3	TN	FA + FC, non-proliferative	y
4	19981001B L	8408	E	YES			1	TN	Focal fibrosis, focal fibrocytic chg	y
5	19981005A R	8408	E	YES			3	TN	microfibroadenoma	y
6	19981008A R	8408	D	YES			3	TN	FC, epithelial hyperplasia, f. atypia	y
7	19981013A L	8408	D	YES			4	TP	DCIS, multifocal, cribif. NG2-S0	y
8	19981013B R	8408	D	YES			4	TN	FC, non-proliferative	y
9	19981014A L	8408	E	YES			4	TN	FC with microcalcifications	y
10	19981022A L	8408	E	NO - NB						
11	19981023A	volunteer	G1	NO - NC						
12	19981023B L	8408	E	YES	TP	TP	5	TP	mxd:DCIS, ID, lobular SIII	y
13	19981028A L	8408	H	YES	TN	TN	4	TN	FC	y
14	19981110A L	post-biopsy	H	NO - NC						
15	19981110B	volunteer	H	NO - NC						
16	19981124A L	8408	H	YES	FN	TP	5	TP	ID and in situ ductal carcinoma	y
17	19981124B L	8408	H	YES	TN	TN	3	TN	FC, FA	y
18	19981125A L	8408	G1	YES	TN	TN	3	TN	FC, proliferative, no atypia	y
19	19981125B R	8408	G1	YES	TN	TN	3	TN	FF, FC, Intraductal papilloma	y
20	19981130A	8408	G1	PENDING					<i>fna only...</i>	
21	19981130B - L-SAR	8408	G1	PENDING					<i>just had biopsy (5-21-99)</i>	
22	19981201A R	8408	G1	YES	TP	TP	4	TP	DCIS+FC (proliferative, +atypia)	y
23	19981202A L	8408	G1	YES	TN	TN	2	TN	FF, apocrine cyst	y
24	19981202A-R	8408	G1	YES	INDET.	FP	2	FP	FF, lobular hyperplasia, adenosis	y
25	19981202B L	8408	G1	YES	FN	FN	5	FN	Invasive ductal +focal DCIS	y
26	19981202B-R	8408	G1	YES	FP	FP	3	FP	FF, epithelial hyperplasia	y
27	19981207A R	8408	G1	YES	TN	TN	3-4A	TN	FC+FAH	y
28	19981207B R	8408	G1	YES	TN	TN	3	TN	Adenosis, fibrosis	y
29	19981209A L	8408	G1	YES	TN	TN	3	TN	FC, FF	y
30	19981209B L	8408	G1	YES	TN	TN	2	TN	FD, Intraductal papilloma	y
31	19981216A L	8408	G1	YES	INDET.	FN	5	TP	ID G2-3 focal tubular features.	y
32	19981221A	8408	G1	YES	TN	TN	3	TN	FC + microcalc.	y
33	19981223A L	8408	G1	YES	TN	TN	3	TN	FC	y
34	19990104A R	8408	G1	NO-WP						y
35	19990104B L	8408	G1	YES	TN	TN	3	TN	FC	y
36	19990106A R	8408	S3	YES	TN	TN	3	TN	FA, FC, FF,	y
37	19990106B R	8408	S2	YES	TN	TN	4A	TN	FF+OSCIFICATION	y
38	19990106C L	8408	G1	YES	TP	TP	5	TP	ID (NOS) poorly differentiated	y
39	19990111A L	8408	G1	YES	TN	TN	3	TN	FC	y
CP	Patient ID	Protocol	Press. P	Acceptable	Blind	Adjunct	BiRads	Unblind	Biopsy	S.conf

40	19990111B R	8408	S3	YES	TN	TN	3	TN	Adenosis		y
41	19990112A L	8408	G1	NO - WP							y
42	19990113A	8408	G1	NO - WP							
43	19990113B- - L-TRE	vol-infl.ca	LS6-RS4	NO - NC							
44	19990120A	8408	Not compl	NO - IP							
45	19990120B L	8408	G1	NO - WP							y
46	19990121A R	8408	G1	YES	TN	TN	4	TN	FC		y
47	19990121B R	8408	G1	NO - WP							y
48	19990121C	8408	LS1-RS1	PENDING					cancelled biopsy only F/U in 6 mths		
49	19990122A L	8408	S4	YES	INDET	TN	3	TN	FC		y
50	19990122A-R	8408	G1	YES	TN	TN	3	TN	FC		y
51	19990125A	8408	LS2-RS2	PENDING							
52	19990127A R	8408		YES	TN	TN	3	TN	Fibroadenoma		y
53	19990127B L	8408	G1	YES	FP	FP	4A	FP	FC+LS in one lobule		y
54	19990128A R	8408	G1	YES	FN	FN	5	FN	FA+DCIS (4mm)		y
55	19990201A L	8408	EUG	NO - EUG							
56	19990203A R	8408	G1	NO - WP							
57	19990203B - H-MAL	8408	EUG	NO - EUG							
58	19990203C R	8408	G1	YES	TN	TN	3	TN	SA		y
59	19990204A R	8408	EUG	NO - EUG							
60	19990204B R	8408	G1	YES	FP	FP	3	FP	FC		y
61	19990208A L	8408	EUG	NO							
62	19990209A L	8408	G1	YES	FN	TP	5	TP	DCIS		y
63	19990209B	8408		PENDING					cancelled biopsy only F/U in 6 mths		
64	19990209C L	8408	S3	YES	INDET	FP	4	FP	FF + PT		y
65	19990209D R	8408	EUG	NO - EUG							
66	19990210A L	8408	EUG	NO - EUG							
67	19990210B-R - I-KAT	8408	G1	YES	TP	TP	4A	TP	ID		y
68	19990210C L	8408	EUG	NO - EUG							y
69	19990216A R	8408	G1	NO - WP							y
70	19990217A R	8408	S2	YES	TN	TN	4	TN	FA (4mm)		y
71	19990218A-R - M-HAL	8408	S6	NO-WP							y
72	19990222A-L - L-SAN	8408	EUG	NO-EUG							
73	19990222B-L - G-BLU	8408	S2	NO-WP-WI							y
74	19990223A-L - A-ROG	8408	S2	YES	FP	FP	3	FP	FF		y
75	19990223B-L	8408	S3	YES	INDET	TP	4	TP	Invasive ductal with tubular features		y
76	19990224A-R - L-ALB	8408	EUG3	NO - EUG							y
77	19990225A-L - M-HAR	8408	E1	NO - WP, WI							y
78	19990301A-R - B-UBI	8408	S2	YES	TN	TN	3,4	TN	Lactacional activity+microcals		y
79	19990301B-L - A-MOR	8408	S2	NO-WP-WI							y
CP	Patient ID	Protocol	Press. P	Acceptable	Blind	Adjunct	BiRads	Unblind	Biopsy	S. conf	
80	19990302A-R - C-COL	8408	G1	YES	FN	FN	5	FN	ID		

81	19990302B R	8408	G1	YES	TN	TN	TN	3	TN FF		y
82	19990302C-L - M-TOR	8408	S4	YES	TN	TN	TN	4A	TN CD+focal atypia+microc+vasculitis sm. Vessels		y
83	19990302D-R - E-WEI	8408	G1	YES	TN	TN	TN	4A	TN FC+ALH+ADH+DH		
84	19990303A-L - G-JOH	8408	S3	YES	INDET	FP	FP	4A	FP ADH+FC(prolif)+microc		y
85	19990303B-L - E-SIL	8408	G1	YES	TN	TN	TN	3	TN Duct dilation with microcal.		y
86	19990309A-R - E-KAS	8408	S3	YES	TN	TN	TN	4A	TN FA		y
87	19990312A - R-PEN			PENDING							
88	19990312B-R - A-PEG	8408	S2	YES	TN	TN	TN	3,4	TN FC+FF+microcal.		y
89	19990316A-R - M-BRU	8408	G1	YES	FP	FP	FP	3	FP FF+FC+AD+DH		y
90	19990316B-R - M-MCG	8408	G1	YES	TN	TN	TN	4A - 4B	TN RS, SA, FA		y
91	19990323A			PENDING							
92	19990324A-R - L-LES	8408	S3	YES	TN	TN	TN	4A	TN CD		y
93	19990325A-L	8408		YES		FN					
94	19990325B-R	8408		YES		TN					
95	19990325C			PENDING							
96	19990326A			PENDING							
97	19990329A-L - G-DAV	8408	805	YES	INDET	NS	NS	4A	NS FC+microcal+sclerosed papilloma		
98	19990330A-R - A-FRA	8408	805	YES	INDET	INDET	INDET	4A	INDET DE+FC+FF+microcal.+FA		y
99	19990330B-R	8408		YES		YN					
100	19990401A-R - B-KIN	8408	G1	YES	FP	FP	FP	4	FP FF+FC		y
101	19990401B-R - A-CHO	8408	G1	NO-MV							
CP											
102	19990405A-R - S-LOV	8408	S3	YES	INDET	TN	TN	3	TN FF+FC		
103	19990406A-R - C-TES	8408	s2	YES	NS	INDET	INDET	5	INDET? DS		
104	19990406B-L - W-FOG	8408	G1	YES	FP	FP	FP	4A	FP DE+FA+microc.		y
105	19990406C-L - A-WHI	8408	S3	NO-WI							y
106	19990407A-R - P-MOU	8408	G1	YES	TN	TN	TN	4B	TN FC non-proliferative		
107	19990412A-R - P-LAF	8408	G1	PENDING					awaiting cc location of lesion		
108	19990412B			PENDING							
109	19990413A			PENDING							
110	19990416A - K-ATH	8408		PENDING							
111	19990419A			PENDING							
112	19990421A-L - M-DIA	8408	S5	YES	TN FP	TN FP	TN FP	4A	FP SA+microc.		
113	19990421B-R - S-COR	8408	G1	YES	FP	FP	FP	4A	FP FC+FA+SA		
114	19990422A-L - E-WEI	8408	G1	YES	TN	TN	TN	4A	TN FN		
115	19990422B-L - R-RIC	8408	S3	YES	TN	TN	TN	US	TN FA		Y
116	19990426A			PENDING							
117	19990426B			PENDING							
118	19990427A-R - A-AAR	8408	S3	YES	NS	FN	FN	5	DCIS+ID		
CP	Patient ID	Protocol	Press. P	Acceptable	Blind	Adjunct	BiRads	Unblind	Biopsy	S. conf	
119	19990427B-R - NO-UTT	8408	G1	YES	TN	TN	4A	TN	FA+FC+LS+microcal		
120	19990428A			PENDING							

COLUMBIA						
Scan Tracking						
ID	EXCLUDED	ACCEPT	PATHOLOGY	AGE	RACE	INTERPRETATION
June-00						
06-01A	11	99	Benign	61	HW	
06-01C	11	99	Benign	35	HB	
06-05A	0	0	Benign	61	HW	TN
06-05B	0	26	Benign	66	HW	
06-07A	0	0	Benign	71	W	TN
06-07B	11	99	Benign	35	HB	
06-07C	26	99	NK		AA	
06-08A	8	99	Benign	52	W	
06-08B	0	0	Benign	50	HW	TN
06-08C	0	27	Benign	64	W	
06-08D	0	30	Benign		W	
06-08E	0	26	Malignant		HW	
06-12A	0	0	Benign	49	W	TN
06-12B	0	0	Benign	56	HW	TN
06-13A L	0	0	Benign	40	HW	TN
06-13B R	0	0	Benign	40	HW	TN
06-14B	0	0	Benign	63	HW	FP
06-15A	26	99	NK		HB	
06-19A	2	99	Benign	39	HW	
06-20A					AA	
06-21A	0	0	Malignant	44	HW	TP
06-22A	0	0	Benign	82	W	TN
06-22B	0	0	Benign	48	HB	FP
06-27A	0	26	Benign	64	W	
06-27B					HW	
06-28A	0	0	Benign	52	HW	FP
06-29A	6	99	Malignant		W	
06-29B	0	0	Benign	50	HW	TN
JUNE TOTALS:						
TOTAL SCANS:						
ACCEPTABLES		(12 benign/1 malignant)				
UNACCEPTABLES						
EXCLUDED PTS						
TBI						

COLUMBIA NEW PP SCANS

July-00					
07-03A	6	99 Benign	46	W	
07-05A	22	99 NK		W	
07-05B	19	99 NK	68	HB	
07-05C	0	17 Benign	59	W	
07-06A L	0	0 Benign	65	HW	FP
07-06A R	0	0 Benign	65	HW	TN
07-11A	0	17 Benign	37	W	
07-11B	11	99 Benign	57	W	
07-13A	0	21 Malignant	54	W	
07-13C	11	99 Benign	38	OTHER	
07-17A	12	99 NK		W	
07-17B	1	99 Benign	43	HW	
07-17C	0	11 Malignant	46	HB	
07-18A L	0	10,11 Benign	49	W	
07-18B R	0	10,11 Benign	49	W	
07-24A	0	11 Benign	38	W	
07-24B			61	W	
07-26A	0	11 Benign	73	OTHER	
07-27B	6	99 Malignant	81	AA	
07-28A	0	0 Benign	57	W	TN
07-31A	6	99 Malignant	66	W	
JULY TOTALS:					
TOTAL SCANS:					
ACCEPTABLES	3 benign				
UNACCEPTABLES					
EXCLUDED PTS					
TBI					

August-00							
08-01A	1	99 Malignant		64	W		
08-01C	11	99 NK		51	W		
08-01D	0	0 Benign		53	W	FP	
08-02A	0	0 Benign		69	HW	TN	
08-02B	0	0 Malignant		46	W	TP	
08-02C	0	0 Benign		44	HB	FP	
08-04A	0	21 Benign		40	HB		
08-09A	26	99 Benign		36	W		
08-09B	1	99 Malignant		76	W		
08-09C	0	0 Benign		54	HW	TN	
08-10A	0	0 Malignant		63	W	TP	
08-10B L	0	13 Benign		55	W		
08-10C R	6	99 Benign		55	W		
08-14A	0	0 Benign		48	W	TN	
08-14B					HB		
08-15A	0	0 Benign		50	AA	TN	
08-15B	0	0 Benign		55	W	FP	
08-16A	11	99 Benign		36	HB		
08-17A	1	99 Malignant		63	W		
08-17B	0	0 Malignant		68	W	TP	
08-21A	0	0 Benign		48	HW	TN	
08-21B	10	99 Atypia		49	ASIAN		
08-21C					W		
08-21D				80	HB		
08-21E	0	21 Benign		60	W		
08-23A	4	99 Malignant		68	HW		
08-28A	11	99 Benign		69	W		
08-29A	3	99 Benign		51	W		
08-29B	0	0 Benign		43	W	TN	
08-29C	0	0 Benign		49	HW	TN	
08-31A	0	0 Benign		61	HW	TN	
08-31B	0	0 Benign		51	W	TN	
08-31C	2	99 Benign		40	HW		
08-31D	0	0 Malignant		82	HB	TP	
08-31E	11	99 Benign		40	HW		
08-31F	6	99 Malignant		54	HW		
08-31G	0	21 Benign		42	HW		
08-31H	4	99 Benign		47	W		
AUGUST TOTALS:							
TOTAL SCANS:							
ACCEPTABLES							
UNACCEPTABLES							
EXCLUDED PTS							
TBI							

COLUMBIA NEW PP SCANS

September-00									
09-05A	4	99	Malignant		80	W	TN		
09-05B	0	0	Benign		52	W	TN		
09-05C	0	0	Benign		60	W			
09-05D	11	99	NK		45	HB			
09-06A	6	99	Benign		79	W			
09-06B	0	0	Benign		43	W	FP		
09-07A	0	0	Malignant		74	HW	TP		
09-07B	0	21	Benign		69	AA			
09-11A	0	0	Benign		54	ASIAN	TN		
09-11B	0	0	Malignant		68	HW	TP		
09-11C	0	0	Benign		70	W	TN		
09-13A						HW			
09-13B						HW			
09-14A	0	0	Malignant		48	W	TP		
09-14B	0	21	Malignant		50	W			
09-14C	7	99	Benign			W			
09-14E	0	0	Benign		46	HW	FP		
09-14F(8)	4	99	Benign		42	HB			
09-18A	0	0	Malignant		80	W	FN		
09-18B	0	0	Benign		42	HB	FP		
09-19A	6	99	Malignant		74	W			
09-19B						HW			
09-20A	0	0	Malignant		67	W	TP		
09-20B						W			
09-20C	0	0	Malignant		75	W	FN		
09-21A	0	0	Benign		63	HB	FP		
09-21B	0	16	Benign		51	W			
09-25A	11	99	Benign		45	W			
09-28A	11	99	Benign		56	W			
09-28B	0	0	Benign		57	HW	FP		
SEPTEMBER TOTALS:									
TOTAL SCANS:									
ACCEPTABLES (9/6)									
UNACCEPTABLES									
EXCLUDED PTS									
TBI									
OVERALL TOTALS									
TOTAL SCANS:									
ACCEPTABLES 117									
UNACCEPTABLES 47(36 benign/11 malignant)									
EXCLUDED PTS 20									
TBI 40									
AVERAGE AGE 10									
55.68									

COLUMBIA NEW PP SCANS

PATIENT SCREENING EXCLUSION CODES

CLINICAL

- 01 Outside of BI-RADS categories 3 and 4
- 02 Not all records available at site for review
- 03 Biopsy at another facility
- 04 Post-menopausal with a palpable lesion
- 05 Pregnant or lactating
- 06 Breast surgery in ipsilateral breast
- 07 Biopsy of ipsilateral breast (core or excisional) within past three months
- 08 Surgical clips or scarring in or on ipsilateral breast
- 09 Patient refused consent.

IMAGING

- 10 Small breasts that, in the judgment of the technologist, cannot be properly positioned
- 11 Lesion is subareolar or located where it cannot be properly illuminated (e.g. lesions visible on MLO view that *cannot also* be visualized on the CC view)
- 12 Patient unable to stand still

CODES FOR UNACCEPTABLE SCANS

BREAST POSITIONING:

- 13 Breast too small
- 14 Soft holder not closed properly
- 15 Non-symmetrical positioning
- 16 Not enough breast tissue in holder

ILLUMINATION

- 17 Size of illumination too small
- 18 Size of illumination too large
- 19 Intensity too high
- 20 Intensity too low
- 21 Inadequate illumination in area of pathology
- 22 Ambient light leakage

DEVICE RELATED

- 23 Computer error
- 24 Air leakage
- 25 Other error: Copy error message number from screen
- 26 Other Device Related: Describe

PATIENT RELATED

- 27 Excessive patient movement
- 28 Patient complaint/refusal after scan has started
- 29 Other Patient Related: Describe

OTHER

- 30 Other: Describe